

# STREAMTASK

## V3.2 Upgrade

There are over 100 corrections and enhancements in StreamTask v3.2. Customers will want to choose a 3-4 week time period between critical EHR updates.

The following is a general timeline for preparation and implementation of the StreamTask version 3.2 upgrade:

- Array and the customer discuss possible timelines for installation, testing and go-live.
- If the current StreamTask instance code has not been updated and/or refreshed to 3.1 in the last 6 months, this should be done. Array reviews disk space and memory on the server and may recommend adjustments including updating IE to v11.
- Array upgrades customer's StreamTask instance to 3.2 and provides testing notes and updated documentation.
- Customer receives System Admin training/review of changes via GoToMeeting.
- Customer reviews and performs testing notes to validate functionality.
- Customer notifies Array that testing is complete.
- If requested, Array provides remote refresher training to testing group as next EHR update begins.