



With the Right Tool, You Can Get the Most Out of Your Updates

StreamTask is the application designed to help manage MEDITECH C/S and 6.X software releases.

StreamTask is a browser-based, platform-independent application that supports software upgrade planning and execution. This new organizational channel allows you to communicate, test, monitor, track and document upgrades with accuracy and confidence.

Organize Your Effort and Streamline Your Process with StreamTask

- Convert DTS Documentation into a database of tasks for testing
- Share an accessible database of testing information and results
- Stay informed through every stage of the update process with real-time status reports
- Enhance communication within your team, and with MEDITECH
- Increase the speed and quality of the testing process
- Track issues throughout their life cycle with detailed documentation
- Meet regulatory requirements with easily available testing documentation

StreamTask enables information sharing, communication, and standardization that **spreadsheets can't provide**. StreamTask is web-based and installed on a local server. StreamTask can also run in a virtual environment.

Organize Your Effort

Download your DTS Documentation from MEDITECH and use built-in conversion tools to import that information into **StreamTask's** SQL database. Now in electronic format, it's easy to quickly share DTSs with end users, IT staff or both. You'll also be able to take advantage of mature task management capabilities that allow you to track, document, communicate and measure your progress.

Streamline the Process

Use **StreamTask's** powerful key word and other search capabilities to analyze the DTS Documentation, and utilize batching tools to find and separate the most significant enhancements from those of less importance, while also eliminating non-applicable tasks in seconds.

Project managers meanwhile can use Browse & Batch to re-assign and re-prioritize workload. Following your hospital's protocols, users document progress and log issues as they occur, thus providing the core team with current, real-time statuses. Users can track Pre-LIVE requirements via Customer Defined Queries and use scripts to facilitate testing procedures. Communication among users is managed within the system, while on-demand and rules-based e-mail messaging keep external staff informed.

About Array Software

Array provides hospitals and health care systems with software tools for successful software upgrade execution and software quality management. Our twenty five years of experience in this field are reflected in our signature upgrade product StreamTask, and in our training and support services.

StreamTask™ Effective Task Management

Request Task Search Reports Dictionary Interface CDR Admin
[Create Save About]

6.0 Core Team, MT 6.06 pp13

Task Type	Completed	Total	Percent
DTSABS	391	521	75
DTSADM	654	691	95
DTSAP	137	143	96
DTSARM	121	121	100
DTSBAR	472	576	82
DTSCA	42	42	100
DTSDR	329	465	71
DTSEDM	296	323	92
DTSEMR	846	1071	79
DTSESS	57	74	77
DTSFA	34	40	85
DTSGL	75	75	100
DTSHIM	504	764	66
DTSHR	397	397	100
DTSITS	492	631	78
DTSLAB	1075	1136	95
DTSMIG	29	35	83
DTSMIS	1185	1481	80
DTSMMA	403	403	100
DTSMRI	198	234	85
DTSMRM	456	481	95
DTSNPR	215	274	78
DTSOM	1639	1776	92
DTSORM	302	499	61
DTSPCM	6	6	100
DTSPCS	2097	2197	95
DTSPHA	836	899	93
DTSPPP	178	215	83
DTSQM	90	138	65
DTSSCH	864	920	94
DTSSS	78	78	100
DTSSYS	111	175	63
DTSUNV	43	66	65
DTSUPT	404	535	76
Totals	15,056	17,482	89

Join us for one of our webcast demonstrations of StreamTask! Go to www.arraysoftware.com/events.html for more information or to register

STREAMTASK The Tool Created to Manage Software Updates

ARRAY

540 Meadow Street Ext., Agawam, MA 01001
Tel: 413-789-2600 ~ Fax: 413-789-2254

On the Web: www.arraysoftware.com